

HIV-1 Reverse Transcriptase, Functional

05-001 200 units, 05-002 1000 units

Shipping and Storage: Ship at 4°C or-20°C and store at -20°C. Do not freezebelow -20°C **Product:** Recombinant full-size HIV-1 Reverse Transcriptase without Tag-peptide expressed in E. coli. It is composed of heterodeimer, p66 and p51 as the one produced in HIV-1 infected cell. It is an RNA-dependent DNA polymerase derived from HIV-1 (AIDS virus), subtype B origin (Ref.1). It also has RNaseH activity and is an enzyme indispensable for reproduction of AIDS virus.

Applications

- 1) It is extremely effective for screening new specific inhibitors for HIV virus as a drug for treating AIDS).
- 2) Generally, Gag and Env proteins are employed as antigens for detecting anti-HIV-1 antibody. However, by using this enzyme in combination as an antigen, the detection will be more sensitive.
- 3) Standards for SDS-PAGE (Fig. 1), Western blotting (Fig. 2), Dot blotting, ELISA

Definition of activity: Activity of intake of 1 nmole of dTMP in 10 min at 37°C is considered as 1 unit using poly(rA) and oligo(dT) as template and primer.

Conditions of measurement: 50 mM Tris-HCl (pH 8.3), 10 mM MgCl₂, 50mM KCl, 3 mM DTT, 0.1% Nonidet P-40, 20 ug/ml poly(rA) · oligo(dT)₁₂₋₁₈, 0.5 mM dTTP ([³H]dTTP,~1 x 10⁵ cpm), and 10-50 units/ml reverse transcriptase.

Activity: ~5,000 units/ml

Purity: Over 90% by SDS-PAGE (CBB staining)

Form: 0.5 mg/ml in 50% glycerol, 40 mM Tris-HCl (pH8.3), 50 mM NaCl, 5 mM MgCl₂, 0.1%

Triton X-100, 1 mM DTT

Data Link: UniProtKB: P04585 GenBank: AAA44988.1




Fig. SDS-Polyacrylamide gel electrophoresis of HIV-1 reverse transcriptase






Fig.2. Western blotting of functional recombinant full-length HIV-1 reverse transcriptase by using anti-HIV-1 Reverse Transcriptase antibody (BioAcademia 65-001).

- 1; 40 ng / lane
- 2; 20 ng / lane
- 3; 4 ng / lane
- 4; 2 ng /lane

Anti-HIV-1 RT antibody was used at 1/2,000 dilution. As second antibody, goat anti-rabbit IgG antibody conjugated with HRP was used at 1/5,000 dilution. ECL system was used.

References: This product was described in Ref 1 and used in the following publications.

- Saitoh A et al. Overproduction of human immunodeficiency virus type I reverse transcriptase in Escherichia coli and purification of the enzyme. Microbiol Immunol 34: 509-521 (1990) PMID: 1699113
- 2.Permanasari ED et al. Enzymatic Activities of RNase H Domains of HIV-1 Reverse Transcriptase with Substrate Binding Domains of Bacterial RNases H1 and H2. Mol Biotechnol. 2015 Jun;57(6):526-38. PMID: 25673083
- 3. Kadokura K et al. Novel urushiols with human immunodeficiency virus type 1 reverse transcriptase inhibitory activity from the leaves of Rhus verniciflua. <u>J Nat Med.</u> 2015 Jan;69(1):148-53. PMID: <u>25349048</u>
- 4.Tada K et al. Abacavir, an anti-HIV-1 drug, targets TDP1-deficient adult T cell leukemia. Sci Adv. 2015 Apr 24;1(3):e1400203. PMID: 26601161
- 5.Izumida M et al. The Spirocyclic Imine from a Marine Benthic Dinoflagellate, Portimine, Is a Potent Anti-Human Immunodeficiency Virus Type 1 Therapeutic Lead Compound <u>Mar Drugs.</u> 2019 Aug 24;17(9). PMID: <u>31450557</u>

Useful Referece: Ref 1 describes infectious cDNA of HIV-1 which was used to construct expression system of this product.

1. Adachi A *et al* "Production of acquired immunodeficiency syndrome-associated retrovirus in human and nonhuman cells transfected with an infectious molecular clone." *J Virol* **59**: 284-291 (1986) PMID: 3016298