

Pfu DNA Polymerase (with dNTPs), Economy

02-021 200 U 2.5U/µl, 02-021-5 200 U 2.5U/µl

Pyrococcus furiosus DNA polymerase (Pfu DNA polymerase) gene was expressed in E. Coli in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and $3' \rightarrow 5'$ exonuclease (proofreading) activity. The MW is 90 kDa, same as that of the natural Pfu DNA polymerase.

■ *Pfu* DNA polymerase is thermostabe and has low error rates.

- It is suitable for PCR and primer extension reactions that require high fidelity synthesis.
- *Pfu* DNA polymerase-generated PCR fragments are blunt-ended.

Applications:

- 1) cloning
- 2) DNA expression
- 3) site-directed mutagenesis

Storage Buffer:

50mM Tris-HCl (pH 8.2), 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.1% Tween20, 0.1% Igepal CA-630,

Store at:

-20°C

- **Concentration:** 2.5 units/ μ l, where one unit is defined as the amount of enzyme that can incorporate 10 nmols of dNTPs into an acid-insoluble material in 30 minutes at 72° C when activated salmon sperm DNA was used as template/primer.
- Quality Assurance: Greater than 95% of protein determined by SDS-PAGE (CBB staining)(Fig.1) The absence of endonucleases and exonucleases was confirmed.
- **PCR Test:** Good amplification result was obtained in PCR reaction using λ DNA as a template (Fig.2).

Reagents Supplied with Enzyme:

- (1) 10 x Standard Buffer (Pfu): 200mM Tris-HCl (pH 8.8), 100mM KCl, 100mM (NH4)2SO4, 20mM MgSO₄, 1% TritonX-100, 1 mg/ml BSA
- (2) 2.5mM (each) dNTPs




Fig.1 SDS-PAGE of Pfu DNA polymerase

Fig.2 Amplification of λ DNA

Related products: # 02-001 Taq DNA Polymerase (+dNTPs) # 02-011 Taq DNA Polymerase BioAcademia, Inc. Tel. 81-6-6877-2335 Fax. 81-6-6877-2336 info@bioacademia.co.jp http://www.bioacademia.co.jp/en/

General composition of PCR reaction mixture (total 50µl)	
Pfu DNA polymerase (2.5 units/	/ul) 0.5 μl
10 x Standard Buffer (<i>Pfu</i>)	5 µl
2.5mM (each) dNTPs	4 µl
Template	<500ng
Primer 1	$0.2\!\sim\!\!1.0\mu\mathrm{M}$ (final conc.)
Primer 2	$0.2\!\sim\!\!1.0\mu\mathrm{M}$ (final conc.)
Sterile distilled water	up to 50µl