

## Anti- *C. perfringens* collagenase antibody, mouse monoclonal (cp-02)

64-050      100 µg

**Shipping and Storage:** Ship at 4C and store at -20C. Do not freeze.

**Immunogen:** Culture supernatant of *Clostridium perfringens*

**Specific Reactivity:** Reacts with collagenases of *Clostridium perfringens* and *C. histolyticum*

### Applications:

1. Western blotting (1/500~1/1,000)
2. ELISA (assay dependent)

This antibody is useful for detecting food-poisoning *Clostridium* strains.


**Background:** *Clostridium perfringens* is one of the major causative agents of food poisoning. *C. perfringens* produces various gelatinolytic enzymes with molecular masses ranging from approximately 120 to approximately 60 kDa. A gelatinolytic enzyme is present in the largest quantity in the culture supernatant, and this enzyme is purified as collagenase. The collagenase of *Clostridium histolyticum* (68 kDa) is the best studied and characterized.

**Isotype:** mouse IgM

**Product:** 1 mg/ml in PBS, 50% glycerol, filter sterilized.


**Purity:** Purified using Ab-Capture for IgM (ProteNova, Japan)

**Data Link:** UniProtKB: [P43153](#) ((COLA\_CLOPE), [Q46085](#) ((COLH\_HATHI))




**Fig.1. Detection of collagenase of *C. perfringens* by Western blotting with monoclonal antibody (MAb cp-02).**

1. Purified collagenase of *C. histolyticum*
  2. Culture supernatant of *C. perfringens*.
- The 68 kDa band in lane 2 is collagenase of *C. perfringens*.  
The primary antibody was used at 1/1,000 dilution.



**Fig.2.** Titration of antibody reactivity of MAb (cp-02) by indirect ELISA using culture medium of *C. perfringens*.

The wells of plate were coated with culture medium of *C. perfringens* (100  $\mu$ l, 1  $\mu$ g/ml). After blocking with 5% skim milk, 100  $\mu$ l of antibody at the indicated dilutions was added to the each well. HRP-conjugated goat anti-mouse IgG, IgM and IgA (100  $\mu$ l, x2000 dilution) was added. Color was developed with orthophenylenediamine as substrate. Optical densities (OD) measured at 490nm.



**Fig.3.** Titration of collagenase in culture medium of *C. perfringens* by indirect ELISA using MAb (cp-02).

ELISA plate is coated with indicated amounts of the culture medium of *C. perfringens* per well. MAb (cp-02) was used at 1/500 dilution. ELISA was performed as in Fig.2.

**Tale 1.** Immunological reactivity of MAb (cp-02) with various food poisoning bacteria

	ELISA	WB
<i>Clostridium perfringens</i> (ATCC13124)	+	68K
<i>Bacillus cereus</i>	-	
<i>Staphylococcus aureus</i>	-	
<i>Campylobacter jejuni</i>	-	
<i>Salmonella Enteritidis</i>	-	
<i>Vibrio parahaemolyticus</i>	-	
<i>Escherichia coli</i> (ETEC)	-	
<i>E. coli</i> 0157:H7 (EHEC)	-	
Purified Collagenase (from <i>C.histoliticum</i> )	+	68K

**Reference:** There has been no publication using this antibody.