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
A. HeLa cells were incubated with or without (2000 kDa dextran w/o SLO) 0.13 $\mu\text{g}/\text{ml}$ SLO on ice for 5 min. After wash with PBS three times, the cells were further with transport buffer containing propidium iodide at 32° C for 5 min. Semi-intact HeLa cells were incubated with 1.5 mg/ml L5178Y cytosol, an ATP regenerating system, GTP, glucose, and 100 $\mu\text{g}/\text{ml}$ fluorescein-dextran of 3, 10, 40, 70, or 2000 kDa at 32°C for 15 min, and then were resealed by treatment with 1 mM CaCl_2 at 32°C for 5 min. After incubation with DMEM supplemented with FCS for 30 min, the cells were observed by confocal microscopy. Since the cells without SLO treatment did not contain the fluorescence of propidium iodide, differential interference contrast (DIC) image was shown. Bar = 10 μm . B. HeLa cells were treated as described in A, were trypsinized, and were subjected to flowcytometry. The histograms of fluorescein fluorescence of dextran with different molecular weight in PI-positive cells were shown.

Data from Kano F et al

A resealed-cell system for analyzing pathogenic intracellular events: perturbation of endocytic pathways under diabetic conditions. PLoS One. 2012;7(8):e44127.doi:10.1371/journal.pone.0044127.

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1. Walev, I. *et al.* "Delivery of proteins into living cells by reversible membrane permeabilization with streptolysin-O." *PNAS* **98**: 3185-3190 (2001) PMID: [11248053](#) **生きた細胞にタンパク質を導入する方法**
2. Palmer, M. "The family of thiol-activated, cholesterol-binding cytolysins." *Toxicon* **39**: 1681-1689 (2001) PMID: [11595631](#)
3. Maeda, Y. *et al.* GPHR is a novel anion channel critical for acidification and functions of the Golgi apparatus. *Nat. Cell Biol.* **10**: 1135-45 (2008) PMID: [18794847](#) **Permeabilization of cells.**
4. Furukawa K. *et al.* Reduction-triggered fluorescent amplification probe for the detection of endogenous RNAs in living human cells. *Bioconjug Chem.* 2009 May 20;20(5):1026-36. PMID: [19374406](#) **Introduction of probes for RNA into permeabilized human HL cells**
5. Thiery J. *et al.* Perforin activates clathrin- and dynamin-dependent endocytosis, which is required for plasma membrane repair and delivery of granzyme B for granzyme-mediated apoptosis. *Blood* **2010** 115:1582-1593. PMID: [20038786](#) **Promotion of endocytosis**
6. Kano F. *et al.* Hydrogen peroxide depletes phosphatidylinositol-3-phosphate from endosomes in a p38 MAPK-dependent manner and perturbs endocytosis. *Biochim Biophys Acta.* 2011 May;1813(5):784-801. PMID: [21277337](#). **Permeabilization of HeLa cells.**
7. Potez S. *et al.* Tailored protection against plasmalemmal injury by annexins with different Ca²⁺ sensitivities. *J Biol Chem.* 2011 May 20;286(20):17982-91. PMID: [21454475](#) **Permeabilization of HEK cells.**
8. Kano F. *et al.* A resealed-cell system for analyzing pathogenic intracellular events: perturbation of endocytic pathways under diabetic conditions. *PLoS One.* 2012;7(8):e44127. PMID: [22952896](#) **Introduction of molecules into HeLa cells.**
9. Imai A. *et al.* MADD/DENN/Rab3GEP functions as a guanine nucleotide exchange factor for Rab27 during granule exocytosis of rat parotid acinar cells. *Arch Biochem Biophys.* 2013 Aug 1;536(1):31-7. PMID: [23702376](#). **Introduction of antibody into cell**
10. Gao N and Lehrman MA. Mannose-6-Phosphate: A Regulator of LLO Destruction. : Inka Brockhausen (ed.), *Glycosyltransferases: Methods and Protocols*, 2013, Springer, Methods in Molecular Biology, vol. 1022, DOI 10.1007/978- Inka Brockhausen (ed.), *Glycosyltransferases*: 1-62703-465-4_20. Link: [springer.com/protocol/10.1007](https://www.springer.com/protocol/10.1007) **Introduction of mannose-6-phosphate into living cells. (The authors specifically recommend BioAcademia streptolysin O for cell permeabilization)**
11. Matsuto M *et al.* Reconstitution of the targeting of Rab6A to the Golgi apparatus in semi-intact HeLa cells: A role of BICD2 in stabilizing Rab6A on Golgi membranes and a concerted role of Rab6A/BICD2 interactions in Golgi-to-ER retrograde transport. *Biochim Biophys Acta.* 2015 Oct;1853(10 Pt A):2592-609. PMID: [25962623](#) **Introduction of protein (Rab6A) into permeabilized HeLa cells.**
12. Yasuga H. *et al.* Logic gate operation by DNA translocation through biological nanopores. *PLoS One.* 2016 Feb 18;11(2):e0149667. PMID: [26890568](#) **Nanopore formation in bilayer lipid membranes.**
13. Ojima K. *et al.* Myosin substitution rate is affected by the amount of cytosolic myosin in cultured muscle cells. *Anim Sci J.* 2017 Nov;88(11):1788-1793. PMID: [28631391](#) **Permeabilization of cells**
14. Kano F. *et al.* Establishment and phenotyping of disease model cells created by cell-resealing technique. *Sci*

[Rep.](#) 2017 Nov 9;7(1):15167. PMID:[29123170](#)11. **Reversible permeabilization of plasma membrane**

15. Watanabe H. et al. Analysis of Pore Formation and Protein Translocation Using Large Biological Nanopores. [Anal Chem.](#) 2017 Nov 7;89(21):11269-11277. PMID:[28980803](#). **Permeabilization for protein translocation**