

## Anti-L-FNR2 (Leaf Ferredoxin NADP Reductase, isoprotein2) antibody, rabbit polyclonal

81-003 100 μg

**Shipping and Storage:** Shipped at 4°C or -20°C and store at -20°C. Do not freeze.

Immunogen: Purified recombinant maize leaf-FNR2 protein (full-size, no-tag attached)

**Form:** 1 mg/ml in PBS- with 50% glycerol. Filter sterilized. No preservative or carrier protein added.

Purity: IgG, affinity-purified with protein A agarose.

Reactivity: Plant L-FNR2 proteins including that of maize and arabidopsis.

Reacts also with other FNR isoforms, Maize L-FNR1 and L-FNR3, and Arabidopsis L-FNR1

## Applications:

- 1. Western blotting (1/2,000-1/50,000 dilution)
- ELISA (assay dependent)
  Other applications have not been tested.

**Background:** Ferredoxin-NADP reductase, leaf isozyme 2 (L-FNR2) plays a key role in regulating the relative amounts of cyclic and non-cyclic electron flow to meet the demands of the plant for ATP and reducing power.

Subcellular location: Chloroplast

Data Link: Swiss-Prot Q8W493 (A. thaliana), Q9SLP5 (Z. mays)




Fig.1 Western Blot of L-FNR2 in plant leaf extract

Anti-L-FNR2 antiserum was used at 1/2,000 dilution. Second antibody (goat anti-rabbit IgG antibody HRP-conjugated, ab97051) was used at 1/10,000 dilution.

- 1. A. thaliana leaf extract (10 µg)
- 2. Z. mays leaf extract (10 µg)

As shown, this antibody cross-reacts with other L-FRN isoproteins. The molecular masses of mature forms of maize FNR1, FNR2 and FNR3 are 34.97, 35.57 and 34.7 kD, respectively.

BioAcademia, Inc. Tel. 81-6-6877-2335 Fax. 81-6-6877-2336 info@bioacademia.co.jp http://www.bioacademia.co.jp/en/



## BSC MC




Fig.2 Cellular distribution of maize FNR isoproteins detected by western blotting

BSC: Proteins (4  $\mu$  g) extracted from bundle sheath cells.

MC: Proteins (4  $\mu$  g) extracted from mesophyll cells.

Anti-FNR2 antibody was used at 1/50,000 dilution

**References**: This product has been used in the following publications.

1. Twachtmann M, Altmann B, Muraki N, Voss I, Okutani S, Kurisu G, Hase T, Hanke GT. "N-terminal structure of maize ferredoxin:NADP+ reductase determines recruitment into different thylakoid membrane complexes. Plant Cell. 2012, Jul;24(7):2979-91. PMID: 22805436 WB; Maize, Arabidopsis