

Anti- Ferredoxin mix (plant) antibody, rabbit polyclonal 81-015 100 µg

Shipping and Storage: Shipped at 4°C or -20°C and store at -20°C. Do not freeze.

Immunogen: A mixture of four Maize Ferredoxin isoproteins, Fd1, Fd2, Fd3 and Fd4.

Form: 2 mg/ml in PBS- with 50% glycerol. Filter sterilized. Azide- and carrier-free.

Purity: IgG, affinity-purified with Protein A.

Reactivity: Essentially all plant Ferredoxin (Fd) isoproteins including those of Arabidopsis and Maize.

Applications:

- 1. Western blotting (1/1,000-1/10,000 dilution)
- 2. ELISA (assay dependent)

Other applications have not been tested.

Background: Ferredoxins are iron-sulfur proteins that transfer electrons in a wide variety of metabolic reactions. Occupies a key position both for transferring the photoreducing power to Fd-NADP+ oxidoreductase (FNR), hence the formation of NADPH, and for mediating the cyclic electron flow around photosystem I (PSI).

Data Link: Swiss-Prot <u>004090</u> (A. thaliana Fd1), <u>P27787</u> (Z. mays Fdx1)




Fig.1 Western Blot of Fd in plant leaf extract.

Anti-Fd antiserum was used at 1/1,000 dilution. Secondary antibody (goat anti-rabbit IgG antibody HRP-conjugated, ab97051) was used at 1/10,000 dilution.

- 1. Arabidopsis leaf extract, 10 µg
- 2. Maize leaf extract, 10 µg

Molecular masses of Fd isoproteins are about 12 kDa, but they migrate at the position around 16-17 kDa on the SDS-PAGE gel due to their strong acidic nature.

1. Hase T. et al. Molecular cloning and differential expression of the maize ferredoxin gene family. Plant Physiol. 1991 May:96(1):77-83. PMID: 16668188 WB; maize

BioAcademia, Inc. Tel. 81-6-6877-2335 Fax. 81-6-6877-2336 info@bioacademia.co.jp http://www.bioacademia.co.jp/en/