

## Anti-Diphtheria Toxin antibody, mouse monoclonal (cdt-01)

Size   100 µg	Product code	64-012
Concentration  Buffer PBS with 50% glycerol  Purity Purified IgG fraction with protein A from hybridoma cell culture medium  Immunogen Toxoid (Diphtheria toxin derived from PWB strain)  Isotype mouse IgG k  Reactivity React with Diphtheria toxin  Validation Specificity has been validated by western blotting  1.Western blotting: x1/400-800 (Fig.1)  Background Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1·193) and B (37 kDa, aa194·535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Size	100 μg
Buffer PBS with 50% glycerol  Purity Purified IgG fraction with protein A from hybridoma cell culture medium  Immunogen Toxoid (Diphtheria toxin derived from PWB strain)  Isotype mouse IgG x  Reactivity React with Diphtheria toxin  Validation Specificity has been validated by western blotting  Application 1.Western blotting: x1/400-800 (Fig.1)  Background Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Storage	-20℃
Purity         Purified IgG fraction with protein A from hybridoma cell culture medium           Immunogen         Toxoid (Diphtheria toxin derived from PWB strain)           Isotype         mouse IgG κ           Reactivity         React with Diphtheria toxin           Validation         Specificity has been validated by western blotting           Application         1. Western blotting: x1/400-800 (Fig.1)           Background         Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Concentration	1.0 mg/ml
Immunogen         Toxoid (Diphtheria toxin derived from PWB strain)           Isotype         mouse IgG κ           Reactivity         React with Diphtheria toxin           Validation         Specificity has been validated by western blotting           Application         1.Western blotting: x1/400-800 (Fig.1)           Background         Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (C-domain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Buffer	PBS- with 50% glycerol
Reactivity   React with Diphtheria toxin	Purity	Purified IgG fraction with protein A from hybridoma cell culture medium
Reactivity React with Diphtheria toxin  Specificity has been validated by western blotting  1.Western blotting: x1/400·800 (Fig.1)  Background  Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypetide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1·193) and B (37 kDa, aa194·535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Immunogen	Toxoid (Diphtheria toxin derived from PWB strain)
Application  1. Western blotting: x1/400·800 (Fig. 1)  Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1·193) and B (37 kDa, aa194·535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Isotype	mouse IgG к
Application  1.Western blotting: x1/400-800 (Fig.1)  Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Reactivity	React with Diphtheria toxin
Background  Corynebacterium diphtheriae infects the nasopharynx or skin. Toxigenic strains secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Validation	Specificity has been validated by western blotting
secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (C-domain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.	Application	1.Western blotting: x1/400-800 (Fig.1)
secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (C-domain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to penetrate the host cell.		
Data Link UniProtKB: Q5PY51 CORDP	Background	secrete a potent exotoxin which may cause diphtheria. The symptoms of diphtheria include pharyngitis, fever, swelling of the neck or area surrounding the skin lesion. Diphtheria toxin is composed of a single polypeptide chain of 535 amino acids (58 kDa) and nicked by cellular protease to give fragments A (Cdomain, 21 kDa, aa 1-193) and B (37 kDa, aa194-535) which are linked by disulfide bridges. Binding to the cell surface of fragment B allows fragment A to
CILL LOVING GOT TOT CONTIN	Data Link	UniProtKB: Q5PY51_CORDP
Please note: All products are FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC		

PROCEDURES. NOT FOR MILITARY USE.



Data Images: 64-012 Anti-Diphtheria toxin antibody, mouse monoclonal (cdt-01)




Fig.1. Western blotting (WB) of cdt-01 antibody.

Diphtheria toxin (BioAcademia, 01-517, 1: 1μg, 2: 0.5μg) was applied to SDS-PAGE and WB. The antibody was used at 1/400 dilution. The HRP-conjugated goat anti-mouse IgG was used at 1/4,000 as the second antibody and visualized by DAB (3,3'-Diaminobenzidine). A protein band 58 kDa in size corresponds to the expected size of the diphtheria toxin.

References: This antibody has not yet been used in publication.

## Related products:

01-517 Diphtheria toxin, functional

64-010 Anti-Diphtheria Toxin, rabbit serum

64-011 Anti-Diphtheria Toxin antibody, rabbit polyclonal